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{-----}

News about deep breathing tests



The determination of the respiratory sinus arrhythmia by means of paced breathing is the most important test of the vagus function. With the BioSign concept of flexibility, tone and dynamics, the individual aspects of vagus regulation can be recorded and assessed in a differentiated manner: A well-functioning vagus should ensure a sufficient basic tone (low heart rate) and be able to adapt over a wide control range (large oscillations of the heart rate = high flexibility) and regulate quickly (rapid drop in heart rate = high dynamics).

There are now indications that the drop in heart rate when breathing out is of particular importance. In a study published in 2016 on 941 heart attack patients, the sinus arrhythmia caused by exhalation proved to be a strong prognostic marker with regard to the survival rate (hazard ratio: 3.4).

Expiration-Triggered Sinus Arrhythmia Predicts Outcome in Survivors of Acute Myocardial Infarction. Sinnecker D, Dommasch M, Steger A, Berkefeld A, Hoppmann P, Müller A, Gebhardt J, Barthel P, Hnatkova K, Huster KM, Laugwitz KL, Malik M, Schmidt G. J Am Coll Cardiol. 2016 May 17;67(19):2213-2220

However, it was not only modern computer-aided medicine that was able to establish the connection between exhalation, HRV and importance for health. In fact, the first known description of the connection between a drop in heart rate during exhalation and health can already be found in the "The Yellow Emperor's Classic of Medicine", the standard work of the traditional Chinese medicine: "Counting four or more pulse beats per exhalation is a sign of death."

In the past 20 years we have seen countless deep breathing tests at BioSign. Our observations support the assumption that cardiac regulation during exhalation is a particularly good indicator of vagal function: the better the physical level of training, the faster the heart rate drop is usually during exhalation.


Heart rate curves during the deep breathing test with paced breathing. Left CHD patient, right competitive athlete. Characteristic is the very rapid and strong drop in heart rate when the athlete exhales compared to the slow and low drop in the CHD patient.

Reason enough to pay more attention to the heart rate regulation when exhaling in the latest version of the HRV scanner, especially since the publications by Guzik and Piskorski on "Heart Rate Asymmetry" provide a suitable framework and, above all, understandable and descriptive parameters.

Heart rate asymmetry by Poincaré plots of RR intervals. Guzik P, Piskorski J, Krauze T, Wykretowicz A, Wysocki H. *Biomed Tech (Berl)*. 2006 Oct;51(4):272-5.

New parameter Guzik index

The heart rate increases when you breathe in and decreases when you breathe out. The RR intervals behave inversely: When the heart rate drops, the RR intervals lengthen. In the Poincaré diagram, inhalation and exhalation can therefore be distinguished: All points above the bisector correspond to a slowing down of the heart rate (increase in RR intervals) and can be assigned to exhalation, all points below the bisector express an acceleration of the heart rate and fall in the inhalation phase.

Poincaré diagram of an deep breathing test of a competitive athlete. up_i indicates the distance between the heartbeat (i) and the bisector during exhalation, $down_j$ the distance between the heartbeat (j) and the bisector during inhalation.

The average distance of all points from the bisector is the well-known SD1. According to Guzik, the SD1 can be represented as the sum of the SD1 values above (exhale) and below (inhale) the bisector:

$$SD1^2 = \frac{1}{n} \sum_{i=1}^n D_i^2 = \frac{1}{n} \sum_{j=1}^{n_{up}} u_j^2 + \frac{1}{n} \sum_{k=1}^{n_{down}} d_k^2$$

$$= SD1_{up}^2 + SD1_{down}^2,$$

If you put the two amounts above and below the bisector in relation to each other, you get the Guzik index:


$$\text{Guzik-Index} = \frac{SD1_{up}^2}{SD1_{down}^2}$$

The Guzik index is the higher, the stronger and faster the heart rate drops during exhalation and, in our view, represents a useful addition to the generalized parameter SD1. With the help of the Guzik index, the dynamic component of vagal regulation can be viewed in a differentiated manner. The Guzik index is also helpful in assessing whether the dynamics are “healthy”: In the case of arrhythmia or “erratic” sinus rhythm (heart rate fragmentation) that is more unhealthy, SD1 and Guzik index change in opposite directions, see table.

Description	Parameter	Ranking (%)	Rank-Chart
Athlete	SD1: 134,58 ms	99,5	
	Guzik-Index: 3,98	94,82	
Atrial fibrillation	SD1: 149,95 ms	99,8	
	Guzik-Index: 1,05	7,2	
Erratic Sinus Rhythm	SD1: 94,45 ms	99,2	
	Guzik-Index: 0,99	5,3	

We have analyzed approx. 800 RSA deep breathing test from our pool for the Guzik index and stored age-corrected reference values so that you can see at a glance whether there are abnormalities in the distribution of the dynamics between inhalation and exhalation.

We have also implemented the parameter in the analysis of the Qiu measurements. It is probably interesting for one or the other user to see whether and how much the dynamics can be trained during exhalation. The following example shows the course of the Guzik index over 12 months. The improvement is significant.

Significantly positive course of the Guzik index over 12 months of Qiu training ($p < 0.01$)

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BioSign GmbH